

## **Title: SARS-CoV-2 genomes recovered by long amplicon tiling multiplex approach using nanopore sequencing and applicable to other sequencing platforms**

**Authors:** Paola Cristina Resende <sup>1,2,5</sup>, Fernando Couto Motta <sup>1</sup>, Sunando Roy <sup>2</sup>, Luciana Appolinario <sup>1</sup>, Allison Fabri <sup>3</sup>, Joilson Xavier <sup>4</sup>, Kathryn Harris <sup>5</sup>, Aline Rocha Matos <sup>1</sup>, Braulia Caetano<sup>1</sup>, Maria Orgeswalska<sup>1</sup>, Milene Miranda<sup>1</sup>, Cristiana Garcia<sup>1</sup>, André Abreu <sup>6</sup>, Rachel Williams <sup>2,7</sup>, Judith Breuer <sup>2, 5, 7\*</sup>, Marilda M Siqueira <sup>1\*</sup>

\* These authors share senior authorship

### ***Filiations***

1 – Laboratory of Respiratory Viruses and Measles, Oswaldo Cruz Institute, FIOCRUZ, Brazil.

2 – University College London, United Kingdom.

3 – Laboratory of Flaviviruses, Oswaldo Cruz Institute, FIOCRUZ, Brazil.

4 – Laboratório de Genética Celular e Molecular, ICB, Universidade Federal de Minas Gerais, Belo Horizonte, Minas Gerais, Brazil

5 – Great Ormond Street Hospital, United Kingdom

6 – National Coordination of Laboratories, Brazilian Ministry of Health, Brasília, Brazil

7 – Pathogen Genomics Unit, University College London, United Kingdom.

### **Summary**

Genomic surveillance has become a useful tool for better understanding virus pathogenicity, origin and spread. Obtaining accurately assembled, complete viral genomes directly from clinical samples is still a challenging. Here, we describe three protocols using a unique primer set designed to recover long reads of SARS-CoV-2 directly from total RNA extracted from clinical samples. This protocol is useful, accessible and adaptable to laboratories with varying resources and access to distinct sequencing methods: Nanopore, Illumina and/or Sanger.

**Keywords:** Coronavirus; Nanopore; GridION; MinION, Illumina, High-throughput sequencing, Long reads, Whole genome, SARS-CoV-2

### ***Short Communication***

The novel Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) belonging to the family of *Coronaviridae* and to the genus *Betacoronaviridae*, emerged in Wuhan, China in December 2019, and has already been introduced in 185 countries to date (1, 2). In UK it has caused more than 166,443 reported cases and 26,097 death and in Brazil 80,246 cases and 5,541 deaths have already been reported, last update April 30<sup>th</sup>, 2020 (2). On March 11th, the WHO declared a SARS-CoV-2 pandemic, reinforcing the need for all countries to implement measures for rapid detection and characterization of the virus to help mitigate virus transmission.

Genomic surveillance has become a useful tool for better understanding virus pathogenicity, origin and spread. Obtaining accurately assembled, complete viral genomes directly from clinical samples is still a challenging task due to the low amount of viral nucleic acid in the clinical specimen compared to host DNA, and to the size of SARS-CoV-2 genome, which is around 30 kb in length. Despite those limitations, we developed a sequencing protocol that successfully obtained whole genomes from SARS-CoV-2 positive samples referred to the National Reference laboratory at FIOCRUZ in Brazil. This protocol was further optimised for higher throughput sequencing at University College London Pathogen Genomics Unit and UCL Genomics to sequence genomes for the COVID-19 Genomics UK Consortium (COG-UK).

The tiling amplicon multiplex PCR method has been previously used for virus sequencing directly from clinical samples to obtain consensus genome sequences (3). This protocol has been applied to Ebola, Zika, Chikungunya and SARS-CoV-2 sequencing (3-6) using preferentially short amplicons (~ 450 base pairs). Nanopore sequencing allows rapid turnaround times (1–2 days) for obtaining a consensus sequence directly from clinical samples and a allows faster response during an outbreak. In this study, we adapted this protocol to recover longer 2 kb reads, decreasing the number of primers required and thus reducing possible mismatches and/or undesired interactions. Additionally, it is easier to assemble larger viral genomes from longer reads enabling higher depth coverage (more than 100 x) in a reduced sequencing time.

Here, we describe three protocols using a primer set designed to sequence SARS-CoV-2 directly from total RNA extracted from clinical samples, which were initially diagnosed using real-time RT-PCR (7, 8). The protocols described herein can be applied to different sequencing platforms, such as Sanger, Illumina and Oxford Nanopore, and therefore are useful, accessible and adaptable to laboratories with different resources and sequencing facilities. By using this protocol, we generated 18 SARS-CoV-2 genomes (15 from clinical samples and 3 isolates from

supernatant of Vero E6 cells culture 72 hours post infection) recovered from different Brazilian states and 22 SARS-CoV-2 genomes from screening healthcare workers in London, UK. The Brazilian genomes are available in the Global Initiative on Sharing All Influenza Data (GISAID) platform.

### ***Primer design***

A set of 17 primer pairs (Supplementary Table 1) was designed to cover approximately 30 kb of the SARS-CoV-2 genome using the online tool primal scheme (<http://primal.zibraproject.org/>) and Prime3 integrated tool in Geneious R9 software, based on the alignment of complete SARS-CoV-2 genomes available in GISAID from the beginning of the outbreak. Each primer pair covers approximately 2 kb of the genome with a 100 bp overlap of amplicons (primer scheme in Table 1).

### ***Primer validation and Sanger sequencing***

The primers were tested in silico using the Geneious R9 software against the 19 available SARS-CoV-2 genomes at that moment. To test the efficiency of each primer pair (10 uM) we performed conventional Sanger sequencing with two positive samples detected in Brazil. cDNA was produced using the SuperScript™ IV First-Strand Synthesis System (Invitrogen) and the DNA was amplified by Q5® High-Fidelity 2X Master Mix (NEB), according to the manufacturer's guidelines. The 17 amplified 2 kb long products were visualized using 1.5 % agarose gel electrophoresis (Figure 1) and each amplicon was sequenced using the BigDye Terminator v3.1 Cycle Sequencing Kit (ThermoFisher, Foster City, USA) and 3.2 μM of the corresponding sequencing primer, indicated in Supplementary table 1. The final reads were recovered by ABI 3130XL Genetic Analyzer (Applied Biosystems) and the sequences were assembled using the Sequencher 5.1 software (Gene Codes).

### ***Genome amplification for Nanopore and Illumina sequencing***

Reverse transcription was initially performed using SuperScript™ IV First-Strand Synthesis System (Invitrogen), using total RNA from samples presenting Ct values  $\leq 30$  for gene E (7, 9). Two multiplexed PCR products (Pool A = 9 primers pairs and Pool B = 8 primers pairs) were generated using the Q5® High-Fidelity DNA Polymerase (NEB) and the primer scheme described in Table 1. The PCR products were purified using Agencourt AMPure XP beads (Beckman Coulter™) and the DNA concentration measured by the Qubit 4 Fluorometer

(Invitrogen) using the Qubit dsDNA HS Assay Kit (Invitrogen). DNA products (Multiplex PCR pools A and B) were normalised and pooled together in a final concentration of 50 fmol.

### ***Library construction for Nanopore sequencing***

The Nanopore library protocol is straightforward as this method is optimised for long reads, such as the generated 2 kb amplicons. Library preparation was conducted using Ligation Sequencing 1D (SQK-LSK109 Oxford Nanopore Technologies (ONT) and Native Barcoding kit 1 to 24 (ONT), according to the manufacturer's instructions. After end repair using the NEBNext® Ultra™ II End Repair/dA-Tailing Module (New England Biolabs, NEB) the native barcodes were attached using a NEBNext® Ultra™ II Ligation Module (NEB). Up to 24 samples were pooled for sequencing in one flow cell. Samples were run in duplicate and two negative controls were used in each round, for experimental validation. A MinION/GridION run using the FLOMIN106 flow cell R9.4.1 was performed comprising 1 to 23 positive samples and 1 negative control (present in all steps post RNA extraction). Sequencing was performed for 2 to 24 hours using the high accuracy base calling in the MinKNOW software. High accuracy base calling was carried out after sequencing from the fast5 files using the Oxford Nanopore Guppy tool. The run was monitored by RAMPART (<https://github.com/artic-network/rampart>) so it could be stopped when  $\geq 20x$  depth for all amplicons was achieved (Figure 2).

### ***Library construction for Illumina sequencing***

As Illumina sequencing chemistry is geared towards sequencing short reads, DNA libraries were generated from the pooled amplicons using Nextera XT DNA Sample Preparation Kit (Illumina, San Diego, CA, USA) according to the manufacturer specifications. The size distribution of the libraries was evaluated using a 2100 Bioanalyzer (Agilent, Santa Clara, USA) and the samples were pair-end sequenced (2 x 300bp) on a MiSeq v3600 cycle (Illumina, San Diego, USA).

### ***Data analysis***

Different data analysis pipelines for Illumina and Oxford Nanopore sequencing were used to extract the consensus files from the raw data.

Demultiplexed Fastq files generated from the Illumina sequencing data were used as an input for the analysis. Reads were trimmed based on quality scores with a cutoff of q30 used to remove low quality regions and adapter sequences were removed. The reads were mapped

to Wuhan Strain MN908947, duplicate reads were removed from the alignment and the consensus sequence called at a threshold of 10X. The entire workflow was carried out in CLC Genomics Workbench software version 20.0.

For the Oxford Nanopore sequencing data, the high accuracy base called fastq files were used as an input for analysis. The pipeline used was an adaptation of the artic-ncov2019 medaka workflow (<https://artic.network/ncov-2019/ncov2019-bioinformatics-sop.html>). We used an earlier version of the workflow which used Porechop to demultiplex the reads. The mapping to the Wuhan reference sequence (MN908947) was done using Minimap2 with Medaka used for error correction. This was all carried out within the artic-ncov2019-medaka conda environment (<https://github.com/artic-network/artic-ncov2019>).

### ***Genomic analysis and phylogenetic reconstruction***

To put the genomes from Brazil and UK generated using this protocol in a global context, SARS-CoV-2 genomes from other countries were recovered from GISAID. Any sequences of length less than 29000 nucleotides, having quality issues on GISAID, or where we detected an unusual frameshifting deletion or insertion relative to the SARS-CoV-2 reference sequence were not included in the phylogenetic reconstruction. To identify similar genomes to the genomes produced and not available yet in GISAID we used the CoV-GLUE website (<http://cov-glue.cvr.gla.ac.uk/#/home>), and we used the NextStrain website (<https://nextstrain.org/ncov/global>), in order to observe the topology of genomes already available in GISAID.

The final curated dataset consisting of 122 SARS-CoV-2 genome sequences was aligned using MAFFT (10). Model testing was carried out using jModelTest (11) and the Generalized Time Reversible plus gamma model was identified as best fitting the dataset. A Maximum Likelihood tree (Figure 3) was generated using RaXML (12) with 1000 bootstraps as branch support. Only values greater than 70 were shown on the tree. The tree was visualized and edited using FigTree v1.4.4 (<http://tree.bio.ed.ac.uk/software/figtree/>).

The genomes recovered using the protocol belongs to seven SARS-CoV-2 lineages (A.2 =1, B.1 = 31, B.1.11 = 2, B.1.2 = 2, B.2.2 = 1, B.3 = 2 and B.6 = 1), according to a classification system previously proposed (<http://virological.org/t/a-dynamic-nomenclature-for-sars-cov-2-to-assist-genomic-epidemiology/458>). To subtype the strains by lineage we used the Pangolin version 1 subtyping tool (<https://github.com/hCoV-2019/pangolin>).

## ***Conclusion***

Here we introduce a versatile sequencing protocol to recover the complete SARS-CoV-2 genome based on reverse transcription plus an overlapping long amplicon multiplex PCR strategy, and associated with pipelines to report the data, and recover the consensus files. The protocol was validated with RNA extracted from some of the first COVID-19 cases detected in Brazil and then optimized and developed for automation at two sequencing facilities at UCL (PGU and UCL Genomics) in London UK. Alternative protocols for Illumina platform, based on an initial amplification of larger fragments (8kb and 10.5kb) produced by one-step RT-PCR with high fidelity enzyme blends were also tested. However, they were prone to producing false mutations, likely due errors during amplification. Based on the fact that SARS-CoV-2 remains conserved, presenting few mutations scattered throughout the genome, the possibility of artificial mutations must be ruled out. We have demonstrated that this overlapping long amplicon multiplex PCR protocol suitable for samples with a wide range of viral loads, generating high coverage throughout the viral genome without artificial indels. It worked well on all four platforms tested (MinION, GridION, Illumina and Sanger) making it suitable for labs with distinct expertise, enabling successful rapid sequencing recovery of the SARS-CoV-2 genome directly from clinical samples.

## ***Ethical statement***

The sequencing workflow optimizations were conducted with the purpose of protocol development and the samples used for this optimization were collected as part of the National Brazilian Surveillance and COG-UK London. We did not use any clinical information or any patient data in this study.


## ***Acknowledgments***

We acknowledge the originators of sequences in GISAID ([www.gisaid.org](http://www.gisaid.org) – acknowledgments in supplementary Table 1); the CGLab of Brazilian Ministry of Health and all the Central Laboratories (LACEN) from Brazilian States; the Oswaldo Cruz Institute, FIOCRUZ; Dr. Rivaldo Cunha coordinator of References Laboratories at FIOCRUZ; the UCL Pathogen Genomics Unit; UCL Genomics; Professor Judith Breuer's research group; the Great Ormond Street Hospital and COG-UK London.

## References


1. Coronaviriae Study Group of the International Committee on Taxonomy of V. The species Severe acute respiratory syndrome-related coronavirus: classifying 2019-nCoV and naming it SARS-CoV-2. *Nat Microbiol.* 2020 Apr;5(4):536-44.
2. Dong E, Du H, Gardner L. An interactive web-based dashboard to track COVID-19 in real time. *Lancet Infect Dis.* 2020 Feb 19.
3. Quick J, Grubaugh ND, Pullan ST, Claro IM, Smith AD, Gangavarapu K, et al. Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples. *Nat Protoc.* 2017 Jun;12(6):1261-76.
4. Quick J, Loman NJ, Duraffour S, Simpson JT, Severi E, Cowley L, et al. Real-time, portable genome sequencing for Ebola surveillance. *Nature.* 2016 Feb 11;530(7589):228-32.
5. Xavier J, Giovanetti M, Fonseca V, Theze J, Graf T, Fabri A, et al. Circulation of chikungunya virus East/Central/South African lineage in Rio de Janeiro, Brazil. *PLoS One.* 2019;14(6):e0217871.
6. Quick J. nCoV-2019 sequencing protocol. 2020 [cited; Available from: [https://www.protocols.io/view/ncov-2019-sequencing-protocol-bbmuik6w?comment\\_id=82272&step=19.21](https://www.protocols.io/view/ncov-2019-sequencing-protocol-bbmuik6w?comment_id=82272&step=19.21)]
7. Grant PRT, Melanie A; Shin, Gee Yen; Nastouli, Eleni; Levett, Lisa J. Extraction-free COVID-19 (SARS-CoV-2) diagnosis by RT-PCR to increase capacity for national testing programmes during a pandemic. 2020.
8. Corman VB, Tobias; Brünink, Sebastian; Drosten, Christian; Landt, Olfert; Koopmans, Marion; Zambon, Maria Diagnostic detection of 2019-nCoV by real-time RT-PCR. Berlin, Jan 17th, 2020; 2020.
9. Corman VM, Landt O, Kaiser M, Molenkamp R, Meijer A, Chu DK, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *Euro Surveill.* 2020 Jan;25(3).
10. Katoh K, Standley DM. MAFFT multiple sequence alignment software version 7: improvements in performance and usability. *Mol Biol Evol.* 2013 Apr;30(4):772-80.
11. Posada D. jModelTest: phylogenetic model averaging. *Mol Biol Evol.* 2008 Jul;25(7):1253-6.
12. Stamatakis A. RAxML-VI-HPC: maximum likelihood-based phylogenetic analyses with thousands of taxa and mixed models. *Bioinformatics.* 2006 Nov 1;22(21):2688-90.

**Figure 1 – Electrophoresis of each ~2kb PCR product amplified from a SARS-CoV-2 sample. Products Pool A and Products Pool B correspond to the two multiplex PCR reactions.**




F = Fragment; L = Low DNA Mass Ladder.

**Figure 2 – RAMPART visualization of the run. (A) coverage depth; (B) fragments lengths and (c) time to achieve 20x 100x and 200x of genome coverage.**





**Figure 3 – Phylogenetic tree of 122 SARS-CoV-2 genomes recovered from GISAID and genomes from Brazil (highlighted in green) and UK (highlighted in blue) produced in this study.**



**Table 1 – Primer scheme to recover 2 kilobases amplicon of SARS-CoV-2 genome.**

Pool A		Pool B	
Oligo name	Oligo sequence (5' to 3')	Oligo name	Oligo sequence (5' to 3')
hCoV_F1_2kb	ACCAACCAACTTTTCGATCTCTTGT	hCoV_F2_2kb	CTGCTCAAAATTCTGTGCGTGT
hCoV_R1_2kb	ACACCACCTGTAATGTAGGCCA	hCoV_R2_2kb	GGTCAGCACCAAAAATACCAGCT
hCoV_F3_2kb	AGCGGACACAATCTTGCTAAACA	hCoV_F4_2kb	TTGTGCACTTATCTTAGCCTACTGT
hCoV_R3_2kb	GGTTGTCTGCTGTTGTCCACAA	hCoV_R4_2kb	TGCCAAAAACCACTCTGCAACT
hCoV_F5_2kb	CACTATTGCAACCTACTGTACTGGT	hCoV_F6_2kb	GTACACTGACTTTGCAACATCAGC
hCoV_R5_2kb	CGTGTGTCAGGGCGTAACTTT	hCoV_R6_2kb	AACGGCAATTCAGTTTGAGCA
hCoV_F7_2kb	TGTACGCTGCTGTTATAAATGGAGA	hCoV_F8_2kb	TGGTACAACATTTACTTATGCATCAGC
hCoV_R7_2kb	TTTGACAGCAGAATTGGCCCTT	hCoV_R8_2kb	TGGGTGGTATGTCTGATCCCAA
hCoV_F9_2kb	CCTTGACCAGGGCTTTAACTGC	hCoV_F10_2kb	AGCAAAATGTTGGACTGAGACTGA
hCoV_R9_2kb	ATCATCTACAAAACAGCCGGCC	hCoV_R10_2kb	CCAAGCAGGGTTACGTGTAAGG
hCoV_F11_2kb	GCTGAAATTGTTGCACTGTGAGT	hCoV_F12_2kb	TGCATTCCACACACCAGCTTTT
hCoV_R11_2kb	AGCACCACCTAAATTGCAACGT	hCoV_R12_2kb	TAACAAAGGCTGTCCACCATGC
hCoV_F13_2kb	ACAAAAGAAAATGACTCTAAAGAGGGTTT	hCoV_F14_2kb	CAGGCTGCGTTATAGCTTGAA
hCoV_R13_2kb	TGTGCTACCGCCTGATAGATT	hCoV_R14_2kb	CATGACAAATGGCAGGAGCAGT
hCoV_F15_2kb	TCAGAGTGTGTAAGTGGACAATCAA	hCoV_F16_2kb	ACGTGAGTCTTGTAAAACCTTCTTTTT
hCoV_R15_2kb	GTACCGTTGGAATCTGCCATGG	hCoV_R16_2kb	ACTGCCAGTTGAATCTGAGGGT
hCoV_F17_2kb	GGAATCATCACAACCTGTAGCTGCA		
hCoV_R17_2kb	TAGGCAGCTCTCCCTAGCATTG		

**Table 2 – Information of SARS-CoV-2 genomes available in this study.**

Virus	Acession number GISAID	Lineage*	CT value**	Average Coverage	Percent Genome 20X	Instrument
hCoV-19/Brazil/RJ-314/2020	EPI_ISL_414045	B.1	19.8	136641	99.60	Oxford Nanopore - 24 hours run
hCoV-19/Brazil/BA-312/2020	EPI_ISL_415105	B.1	23.8	137747	99.60	Oxford Nanopore - 24 hours run
hCoV-19/Brazil/ES-225/2020	EPI_ISL_415128	B.2.1	27.7	15914	99.60	Oxford Nanopore - 24 hours run
hCoV-19/Brazil/AL-837/2020	EPI_ISL_427292	B.1	21.0	4524.51	99.60	Illumina MiSeq
hCoV-19/Brazil/BA-510/2020	EPI_ISL_427293	B.1	23.3	6987.78	99.60	Illumina MiSeq
hCoV-19/Brazil/DF-615i/2020	EPI_ISL_427294	B.1	14.4	3626.95	99.60	Illumina MiSeq
hCoV-19/Brazil/DF-619i/2020	EPI_ISL_427295	B.1	13.6	2000.83	99.60	Illumina MiSeq
hCoV-19/Brazil/DF-861/2020	EPI_ISL_427296	B.1	20.3	10283.66	99.60	Illumina MiSeq
hCoV-19/Brazil/DF-862/2020	EPI_ISL_427297	B.1	26.2	3287.77	99.60	Illumina MiSeq
hCoV-19/Brazil/DF-891/2020	EPI_ISL_427298	B.1	24.7	4212.5	99.60	Illumina MiSeq
hCoV-19/Brazil/RJ-352/2020	EPI_ISL_427299	A.2	27.5	3461.93	99.60	Illumina MiSeq
hCoV-19/Brazil/RJ-477/2020	EPI_ISL_427300	B.1	23.2	2526.92	99.60	Illumina MiSeq
hCoV-19/Brazil/RJ-477i/2020	EPI_ISL_427301	B.1	14.6	4557.41	99.60	Illumina MiSeq
hCoV-19/Brazil/RJ-763/2020	EPI_ISL_427302	B.1	19.0	2777.22	99.60	Illumina MiSeq
hCoV-19/Brazil/RJ-818/2020	EPI_ISL_427303	B.1	20.1	5551.05	99.60	Illumina MiSeq
hCoV-19/Brazil/RJ-872/2020	EPI_ISL_427304	B.1	19.9	7583.52	99.60	Illumina MiSeq
hCoV-19/Brazil/SC-766/2020	EPI_ISL_427305	B.6	20.0	6120.73	99.60	Illumina MiSeq
hCoV-19/Brazil/SC-769/2020	EPI_ISL_427306	B.1	19.0	557865	99.60	Illumina MiSeq
hCoV-19/England/LOND-D37F2/2020	not available	B.1	30.0	305.65	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D382C/2020	not available	B.1	26.6	274.079	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D384A/2020	not available	B.1	26.9	232.13	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3868/2020	not available	B.1	27.1	191.063	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3886/2020	not available	B.1	27.2	257.627	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3895/2020	not available	B.1	27.3	189.683	99.51	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D38A4/2020	not available	B.1	28.3	293.545	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D38B3/2020	not available	B.1	28.4	269.598	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D38C2/2020	not available	B.1	28.7	252.471	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D38D1/2020	not available	B.1	25.7	163.62	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D38E0/2020	not available	B.1	24.8	229.004	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D390B/2020	not available	B.1	26.3	394.739	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D391A/2020	not available	B.1	27.2	257.287	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3938/2020	not available	B.1	28.1	298.876	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3947/2020	not available	B.1	28.1	333.132	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3956/2020	not available	B.1	28.7	250.45	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D380E/2020	not available	B.1.11	22.2	325.745	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3929/2020	not available	B.1.11	28.0	587.767	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3877/2020	not available	B.2.1	27.2	262.921	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D38FF/2020	not available	B.2.2	25.7	246.488	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D383B/2020	not available	B.3	26.8	389.535	99.52	Oxford Nanopore - 6 hours run
hCoV-19/England/LOND-D3859/2020	not available	B.3	26.9	333.778	99.52	Oxford Nanopore - 6 hours run

\* Lineage based on Pangolin version 1 subtyping tool (<https://github.com/hCoV-2019/pangolin>); \*\* CT = Cycle threshold , samples from Brazil and UK had the ct value measured by different RT-PCR protocols.