



Uncropped versions of SDS-PAGE gels


PLpro WT Ub  
chain specificity  
assays




S1 site WT/Mts +  
K48 triUb kinetics  
comparison





S1 site WT/Mts +  
proISG15 kinetics  
comparison





PLpro WT +  
K48 triUb  
kinetics assays





PLpro WT +  
proISG15  
kinetics  
assays





PLpro R166S/E167R  
+ K48 triUb kinetics  
assays




PLpro R166S/E167R +  
proISG15 kinetics assays





PLpro N156E +  
K48 triUb kinetics  
assays




PLpro N156E +  
proISG15 kinetics  
assays




PLpro Y171R +  
K48 triUb kinetics  
assays





PLpro Y171R +  
proISG15 kinetics  
assays





PLpro K232E +  
K48 triUb kinetics  
assay





PLpro K232E +  
proISG15 kinetics  
assays



PLpro T225Y +  
K48 triUb kinetics  
assays





PLpro T225Y +  
proISG15 kinetics  
assays





PLpro F69S +  
K48 triUb kinetics  
assays



PLpro F69S +  
proISG15 kinetics  
assays



PLpro WT + ISG15<sup>CTD</sup>  
kinetics assays



PLpro F69S +  
proISG15<sup>CTD</sup> kinetics  
assays

