



# Appendix 3: IBL protocol for setting up the behavioral training rig

The latest version of this protocol can be found [on Figshare via this specific link](#).

All IBL protocols cited in the article can be found [on Figshare via this master link](#).

The latest version of the CAD models and technical drawings for the manufactured rig parts can be found [on Figshare via this specific link](#).

The latest version of the whole rig CAD assembly can be found [on Figshare via this specific link](#).

**Please cite the associated article** when using any of these materials and protocols:

The International Brain Laboratory et al. (2020) *Standardized and reproducible measurement of decision-making in mice*. bioRxiv, 909838. <https://doi.org/10.1101/2020.01.17.909838>