

TOM22/COX2

<i>Number of tests</i>	<i>Calculation time</i>	<i>Colocalization (%)</i>
500 tests	15s	31,66
1000 tests	26s	31,43
2500 tests	58s	31,93
5000 tests	1min56	32,13
1000 tests	4min8s	32,41

TOM22/PMPCB

<i>Number of tests</i>	<i>Calculation time</i>	<i>Colocalization (%)</i>
500 tests	13s	71,29
1000 tests	24s	70,01
2500 tests	59s	69,96
5000 tests	1min50	69,7
10000 tests	3min37	69,11

AURKA/TOMM20

<i>Number of tests</i>	<i>Calculation time</i>	<i>Colocalization (%)</i>
500 tests	13s	24,6
1000 tests	24s	24,8
2500 tests	59s	23,32
5000 tests	1min58s	23,83
10000 tests	3min42s	24,72

AURKA/PMPCB

<i>Number of tests</i>	<i>Calculation time</i>	<i>Colocalization (%)</i>
500 tests	12s	53,31
1000 test	23s	55,13
2500 tests	56s	53,29
5000 tests	1min55s	53,59
10000 tests	3min33s	53,3

AURKA/TOMM20 antigen retrieval

<i>Number of tests</i>	<i>Calculation time</i>	<i>Colocalization (%)</i>
500 tests	21s	19,4
1000 tests	33s	19,74
2500 tests	1min6s	20,47
5000 tests	2min5s	20,86
10000 tests	4min4s	20,4

AURKA/PMPCB antigen retrieval

<i>Number of tests</i>	<i>Calculation time</i>	<i>Colocalization (%)</i>
500 tests	18s	14,63
1000 tests	29s	12,79
2500 tests	1min9s	13,09
5000 tests	1min57s	13,65
10000 tests	3min42s	14,07

